


A Treatment for Everyone and Every Condition

GentleMax Pro[®]



Photos: Kathleen P. Hutton, M.D.



Science. Results. Trust.

Syneron[®]



CANDELA[®]

GentleMax Pro[®]


A Treatment for Everyone and Every Condition

How does it work?

- GentleMax Pro uses laser technology to deliver an intense but gentle burst of energy onto targeted areas of the skin, delivering a controlled amount of targeted therapeutic heat safely and effectively.


How long does treatment take?

- A treatment could take as little as a few minutes to as long as an hour, depending on the area and the condition being treated.




Proven safe and effective for most all skin types.

Laser Hair Removal



Facial and Leg Veins

Sun and Age Spots




What can I expect during treatment?


- Treatment will take place in your practitioner's office. Your practitioner will utilize a small handpiece to deliver laser pulses onto the desired treatment area.
- Some patients may feel a slight discomfort with each laser pulse, sometimes described as a snap of a rubber band on the skin.
- Discomfort is resolved quickly and minimized with the patented Dynamic Cooling Device™ (DCD™) or air cooling. For some patients, a topical anesthetic may be used prior to treatment.

When can I expect to see results?

- The number of sessions depends on the desired treatment and the condition of your skin.




Hair Anatomy





Bulb/bulge are critical structures responsible for hair regrowth.


Melanin Absorption




Wavelength




Laser Hair Removal



Anagen



Catagen




Telogen



Before

Post 3 treatments
Photos: M. Ercan, M.D.

After





Before



After

Post 6 treatments
Photos: O. Zerpa, M.D.



Before



After

Photos: M. Patel






Before



After

Post 5 treatments
Photos: M.C. Kutun, CME, CLS, CLP



Before



After

Post 3 treatments
Photos: M. Elder, M.D.





Sun and Age Spots


What causes them?

- Pigmented lesions are often caused by an abnormal production of melanin, which becomes visible on the skin's surface

How does GentleMax Pro work?


- Long pulsed laser energy delivered at 3 ms closely matches the Thermal Relaxation Time (TRT) of melanocytes
- Heating of the melanocytes in benign pigmented lesions causes thermal damage, resulting in removal of these types of lesions






Facial and Leg Veins

- Sclerotherapy is the *gold standard* for the treatment of leg veins
- The 1064 nm wavelength can clean up what sclerotherapy does not
- Veins less than 2 mm in size respond best



Spider Veins







Syneron and Candela—two extraordinary brands, one global leader

Syneron has a remarkable ability to innovate the newest technologies in aesthetic medicine.
Candela brings four decades of industry-leading scientific expertise.
Together, these two brands offer “gold standard” technologies for advanced aesthetic solutions.

Schedule a personal consultation with your practitioner today!
www.syneron.com | www.candelalaser.com



Science. Results. Trust.

Syneron[®]

 **CANDELA**[®]

© 2012. All rights reserved. Syneron and the Syneron logo are registered trademarks of Syneron Medical, Ltd. Candela, the Candela logo and GentleMax Pro are registered trademarks of Candela Corporation. Dynamic Cooling Device and DCD are trademarks of Candela Corporation and may be registered in certain jurisdictions. PB76861EN-NA